

All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Biomet Inc. or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet. Confidential information intended solely for Zimmer Biomet employees, the Zimmer Biomet sales force and authorized representatives. Distribution to any other recipient is prohibited. For complete product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the package insert and [www.zimmerbiomet.com](http://www.zimmerbiomet.com). Check for country product clearances and reference product-specific instructions for use.

Nothing in this presentation should be construed as a substitute for reviewing the actual study itself. Any reference to the content within this presentation must be accompanied by provision of the actual study.

# Clinical Trials

## Restore and The SternaLock Blu Study

R060518

---

# The Restore Study

**Title:** Sternal Closure with Rigid Plate Fixation Versus Wire Closure: A Randomized Controlled Multicenter Trial

**Objective:** Evaluate sternal bone healing and pain in patients with multiple risk factors undergoing open heart surgery via a median sternotomy.

**Design:** Prospective, randomized, controlled multi-center trial

**Investigational Device:** SternaLock Silver

**Configurations:** SternaLock: 2 X plates & one L plate at the manubrium; Wires: 1 wire per 10 kg and min of 7 wires peristernal.

**Patient Demographics:** High Risk Patients

**Sites:** 6 (USA & Germany)

**Patients:** 140

**Times for follow-up:** Pre-op, during hospital stay, 3 & 6 weeks, 3 & 6 months.

**Primary Endpoints:** Radiographic evaluation of bone healing using CT scans

**Secondary Endpoints:** Pain (visual analog scale), pain medication (index admission)

**Key Results:**

- Improved Healing
  - SternaLock Silver had superior healing at both 3 and 6 months
  - At 3 months SternaLock Silver was 15% sternal union vs. 0% with wire
  - At 6 months SternaLock Silver was 70% sternal union vs. 24% with wires
- Reduced early postoperative pain
  - 25% reduction in pain scores and narcotic usage
- Complication rates
  - Similar between groups



## Data Details:

- Improved Healing
  - Limited data because the same patients were not scanned at 3 & 6 months.
  - Higher risk patients which could impact healing and lead to lower scores when compared to the SLB Study.
- Similar Complication Rates
  - Majority of the StenaLock Silver's complications occurred due to improper screw length (10/12mm screws) used
  - Patients had multiple co-morbidities and high risk for complications
- Pain
  - No correlation between sternal healing and pain due to not scanning the same patient at 3 & 6 months

## Surgeons & Institutions

Restore Surgeons	Restore Institutions	City / State
Jaishankar Raman, MD, PhD	Rush University Medical Center	Chicago, IL
Sven Lehmann, MD	University of Leipzig Heart Center	Leipzig, Germany
Kenton Zehr, MD	Scott and White Clinic	Temple, TX
Brian J. DeGuzman, MD	St. Joseph's Hospital and Medical Center	Phoenix, AZ
Lishan Aklog, MD	St. Joseph's Hospital and Medical Center	Phoenix, AZ
Edward Garrett, MD	Cardiovascular Surgery Clinic, Memphis	Memphis, TN
Michael Wong, MD	University of California Davis Medical Center	Sacramento, CA

### Resource:

Raman J, Lehmann S, Zehr K, et al. Sternal closure with rigid plate fixation versus wire closure: A randomized controlled multicenter trial. Ann Thoracic Surg. 2012; 94:1854-61.

# The SternaLock Blu Study

**Title:** An Evaluation of Rigid Plate Fixation in Supporting Bone Healing: A Prospective, Multi-Center Trial: Final Study Report

**Objective:** To evaluate clinical and economic outcomes following sternal closure with either SternaLock Blu or wire cerclage in patients undergoing elective cardiac surgery.

**Design:** Prospective, randomized, controlled multi-center trial

**Investigational Device:** SternaLock Blu

**Configurations:** SternaLock Blu: 2 X plates & one (L, JL, hexagon or box) plate at the manubrium; Wires: min of 6 wires in any configuration per institutional preference

**Patient Demographics:** Elective Cardiac Surgery Patients

**Sites:** 12 (USA)

**Patients:** 236

**Times for follow-up:** Pre-op, during hospital stay, 3 & 6 weeks, 3 & 6 months.

**Primary Endpoints:** Radiographic evaluation of bone healing using CT scans


**Secondary Endpoints:** Pain (visual analog scale), pain medication (thru 6 months), Quality of Life (SF-36), Function Outcomes (UEFI), Complications, Health Economics.

## Key Results:



- **Improved Healing**

- Sternal healing was assessed with CT scans by independent radiologists using validated methods.
- At 3 months SLB was 2.6X higher vs. wires
- At 6 months SLB was 1.2X higher vs. wires



- **Fewer Sternal Complications**

- 0% complication rate SLB vs. 5% complication rate wires
- Reflects all sternal complications over 6 months (vs. "<1%" rate often cited as the 30 day DSWI rate)
- In a regression analysis, the method of closure was the only predictor of sternal complications.



- **Improved Patient Recovery**

- 237 total fewer days in rehabilitation and recovery for SternaLock Blu patients compared to wire patients. \*

**Proven Cost Savings**



- **Cost Effective Solution**

- Proven cost savings over 6 months with SternaLock Blu.
- Increased costs with wire patients primarily due to increased sternal complications and greater outpatient resource utilization
  - (e.g. rehab hospitals, skilled nursing facilities).



**Data Details:**

- **Improved Healing**
  - Higher rates of healing and faster healing with SLB
- **Sternal Complications**
  - Probability of complications with SLB remains low even with increasing BMI
  - Sternal closure was the only factor that was found to be a predictor for complications and infections
  - Patients with sternal complications (n=6) required 11 reoperations and 94 days of additional hospital stays
- **Proven Cost Savings**
  - Primary drivers of cost savings were cost related to sternal complications and outpatient resource utilization (OPRU) (eg. 237 fewer days)

**Surgeons & Institutions**

SLB Study Surgeons	SLB Study Locations	City/State
Keith B. Allen	St. Luke's Mid America Heart Institute	Kansas City, MO
Vinod H. Thourani	Emory University	Atlanta, GA
Yoshifumi Naka	Columbia University Medical Center	New York, NY
Kendra J. Grubb	University of Louisville	Louisville, KY
John Grehan	United Heart and Vascular Clinic	Saint Paul, MN
Nirav Patel	Lenox Hill Hospital, New York, NY	New York, NY
T. Sloane Guy	Temple University, Philadelphia, PA	Philadelphia, PA
Kevin Landolfo	Mayo Clinic, Jacksonville, FL	Jacksonville, FL
Marc Gerdisch	Franciscan St. Francis Health, Indianapolis, IN	Indianapolis, IN
Mark Bonnell	University of Toledo, Toledo, OH	Toledo, OH
David J. Cohen	St. Luke's Mid America Heart Institute	Kansas City, MO
Donald Botta	Florida Hospital	Orlando, FL

**Resources**

This is a summary of the Study entitled, CR 0712: An Evaluation of Rigid Plate Fixation in Supporting Bone Healing: A Prospective, Multi-Center Trial; Final Study Report CR0712S & CR0712E. Report on file at Zimmer Biomet.

\*SternaLock Blu patients spent a total of 237 total fewer days in rehab hospitals or skilled nursing facilities over 6 months (705 total days for wire cerclage (n=120) vs. 468 total days for SternaLock Blu patients (n=116).